



TECHNICAL DATA SHEET

Ami Polymer Pvt. Ltd.
"Sealing Expert in Silicone"

APPROVED DOCUMENT

Ami Polymer Pvt Ltd clean room molded PP Pinch clamp is designed to address the needs of the single use disposable and biopharmaceutical industries. Pinch clamp used to control the fluid flow rate. We produce pinch clamps to be used with flexible tubes that work in low-pressure application. Pinch Clamp manufacturing process carried out in clean room of ISO Class 7 where the process is cautiously controlled from the raw material to the finished product manufactured.

FEATURES/BENEFITS

- Biopharmaceutical grade
- Easy installation
- Autoclavable and Gamma sterilisable
- MOC: Polypropylene

APPLICATION

- Single Use System
- Vaccine and Bio therapeutic production
- Fluid flow control



CERTIFICATION

- FDA 21 CFR 177.2600
- USP Class VI

ORDERING INFORMATION

Part No.	Size	Tube OD (in inch)
PP-PC-0250S	Small	1/8" - 1/4"
PP-PC-0375M	Medium	1/4" - 3/8"
PP-PC-0500L	Large	1/2" - 3/4"



For, Ami Polymer Pvt. Ltd.
Quality Assurance

"ELECTRONICALLY APPROVED DOCUMENT"

This information is believed to be accurate and refers to the laws, regulations and products at the date of issue. However It is the responsibility of our customers to determine that their use of Ami Polymer Pvt. Ltd.'s products is safe, lawful, and technically suitable for their applications. Because of possible changes in the laws and regulations, we cannot guarantee that the status of this product will remain unchanged.